

FOCUS

DENTAL UNIT



D E S I G N & F U N C T I O N I N H A R M O N Y

Qdent

FOCUS

DENTAL UNIT



Danish Design at its best

Danish design is known for its ability to combine optimum functionality and visual delight. With the Focus unit, Qdent will continue this tradition for the daily benefit of the users. We believe that a beautiful, timeless design is an important factor. After all, a dental unit is a considerable investment and you should enjoy looking at it for many years.

Made in the Netherlands

The focus Unit is nowadays manufactured and further developed in the Netherlands. The Netherlands, being a small country, have always depended on inventiveness and their high standard of technology, for their economic survival. With the production of the Focus Unit this will not be any different.

QUALITY AND FLEXIBILITY We are never satisfied until our customers are. Our philosophy with regard to design, development and manufacturing of user-friendly dental equipment is really that simple. In order to ensure that we always meet the needs and expectations of our users, we use nothing but *the best materials*. We develop the Focus concept on a continuous basis in close collaboration with our customers and distributors and we ensure a *flexible response* at all times to customers' special requirements through comprehensive proprietary production.

DEVELOPMENT AND RELIABILITY All new components are subjected to stress tests far in excess of normal dentistry equipment standards. They are then tested under everyday conditions by our clinical test teams. Behind the scenes, elegant component design, strong materials and straightforward electronic control mean long-lasting operational reliability, while the leading role is played by the Focus concept of a functional and reliable workplace that meets all *the advanced standards* of today with regard to quality, working comfort and functionality.

FOCUS



1 BEYOND ALL EXPECTATIONS

WITH FOCUS 1 you are always on top. Form, design and details create a restful atmosphere and bring functional harmony. Everything is within reach, including as many as 6 instruments – each suspended in its own spring-balanced holder to prevent unnecessary muscle strain.

The *ergonomic design* of the chair ensures a relaxed working posture and its soft lines make it comfortable and relaxing for the patient. Focus can easily be adjusted during treatment. The versatile instrument panel means that instruments can always be placed just where you need them. You can program the electronic control system to ensure, for example, that the chair is automatically adjusted to the position you prefer. The upper surface of the instrument panel features a number of touch-sensitive foil contacts for easy operation of fountain, cup filler, bowl flush, assistant call and automatic cut in/out of the chip blow.

FOCUS 1 is designed to be *attractive and functional*, available in various exciting colours.



FOCUS

2 DESIGN AT WORK



FOCUS 2 merges *quality and functionality in an attractive design* – smooth adjustment and use, a light touch to move the instrument panel and programmable electronic control to make sure the chair is always adjusted correctly. A built-in safety function ensures that the movement of the chair ceases automatically when the seat or back encounter resistance. The back is so thin that the clinical assistant or dentist can easily place his or her thigh and knee under the back in order to attain a *comfortable working position*. For the patient, the Focus chair is comfortable and relaxing. The concave design of the headrest with its raised, rounded edges offers firm and secure support for the patient's head. Both armrests can be swung to the side, providing unrestricted access for the patient – or you can even remove them completely.

FOCUS 2 has room for up to 6 suspended instruments and comes with a practical, corrosion-free tray. The unit features touch-sensitive foil contacts on the front of the instrument panel.

FOCUS 2 is a functional and comfortable workplace and a source of colour for every surgery.



FOCUS NUMEROUS DETAILS

The Focus offers numerous details in both design and functionality. Just some are published in this brochure. Please contact your local dealer for exact specifications of the unit and a list of options & possibilities.



Easy-to-clean spittoon bowl

The hygienic spittoon has no sharp edges or joints and is easy to clean and dismount.



The unique ball swivel joint system

This system is the epitome of design simplicity. Yet it is more. The recently developed ball swivel joint system ensures optimum resistance-free flexibility of the instrument arm.



Focus aspiration holder

The Focus aspiration holder comes either suspended from the instrument panel or mounted on the unit column. The holder mounted on the unit column has a telescopic function and easily slides towards the patient. You can also choose for a combination of the two options. The aspiration adaptors are easily removed for cleaning and disinfection.



Instrument tray suspension with 1 or 2 trays

The instrument tray is easily positioned to the left or right of the patient. After treatment it is easy to swing the tray back. Choice of 1 or 2 trays.

WORKING COMFORT & DESIGN



A long day at the clinic makes freedom of movement and operation an absolute necessity. Qdent has focused on every detail of every function to ensure that instrument holders can be positioned in the best place, depending on the *preferred work posture*. Even the chair design enables the dentist to give his or her best. The soft comfort of the chair also makes things as pleasant as possible for the patient.

The soft lines and rounded shapes of the Focus embrace modern trends with a design that is appropriate for all treatment situations. The *combination of technology and design* also reassure the patient that he or she are in safe hands when they make themselves comfortable in the chair.



CHOOSE YOUR COLOR

Kelvinring 5
2952 BG Alblasterdam
The Netherlands


t +31 (0)78 68 40 954

f +31 (0)78 68 49 649

i www.qdent.nl

e info@qdent.nl

distributed by:



The logo for Qdent, featuring a stylized yellow 'Q' followed by the word 'dent' in a grey serif font.

